

## **Influenza A and B virus HA/NA/(A-MP) amplification**

This protocol is for HA, NA and MP amplification for influenza A, or for HA and NA genes amplification of influenza B. PCR products for both influenza A and B can be amplified in one reaction using FluA or FluB primers and SuperScriptIV One-Step RT-PCR System (ThermoFisher Cat#12594100). The PCR amplicons are then quantified for further library preparation for NGS using chosen library preparation method for NGS.

### **Primers:**

Target	Primer Name	Primer sequence	Pooling Vol (µl)
FluB- HA/NA	B-HANA-UniF	GGGGGGAGCAGAAGCAGAGC	13
	B-HANA-UniR	CCGGGTATTAGTAGTAACAAGAGC	13
FluA- HA/NA/MP	A-HA-UniF	GGGGGGAGCAAAAGCAGGGGA	5
	A-HA-UniR	CCGGGTATTAGTAGAAACAAGGGTG	5
	A-NA-UniF	GGGGGGAGCAAAAGCAGGAGT	6
	A-NA-UniR	CCGGGTATTAGTAGAAACAAGGAGT	6
	A-M-UniF	GGGGGGAGCAAAAGCAGGTAG	2.5
	A-M-UniR	CCGGGTATTAGTAGAAACAAGGTAG	2.5
Total pool (50 reactions)			53

### **RT-PCR using SSIV RT-PCR kit**

#### **1. Make Master Mix in a clean room.**

Master Mix	Volume (µl)	
	x 1	x 10
Primer Pool	1	10
2 x Master Mix	12.5	125
Reverse Transcriptase	0.25	2.5
Water	3.25	32.5
<b>Total</b>	<b>17</b>	<b>170</b>

Label a 1.5 ml tube, make master mix by adding reagents in the above table in the tube. Add Reverse Transcriptase at the last step and vortex briefly to mix and then a quick spin to collect all contents at the bottom of the tube.

## 2. Add RNA to Master Mix

Take the tube into the extraction room, aliquot **17 µl** of master mix into each PCR tube.

Then add **8µl** of RNA to each PCR tubes containing master mix, change tips for every well. The total volume should be **25µl**.

A quick vortex to mix and then spin the tubes briefly to collect all contents at the bottom of the tube.

Put in the PCR machine to start the RT-PCR program:

RT-PCR	
FluA/FluB NGS	
SSIV RT-PCR	
55°C	10 minutes
98°C	2 minutes
98°C	10 seconds
45°C	10 seconds
72°C	2 minutes
x5	
98°C	10 seconds
57°C	10 seconds
72°C	2 minutes
x35	
72°C	5 minutes
4°C	∞